

# *Self-Supervised Deep Learning Reconstruction for Highly Accelerated Diffusion Imaging*

*Ismail Arda Vurankaya<sup>1</sup>, Yohan Jun<sup>2,3</sup>, Jaejin Cho<sup>2,3</sup>, Berkin Bilgic<sup>2,3,4</sup>*


*<sup>1</sup> Bogazici University, Istanbul, Turkey*

*<sup>2</sup> Athinoula A. Martinos Center for Biomedical Imaging, Charlestown, MA, United States*

*<sup>3</sup> Department of Radiology, Harvard Medical School, Boston, MA, United States*

*<sup>4</sup> Harvard/MIT Health Sciences and Technology, Cambridge, MA, United States*

# Motivation






Multishot EPI

- Multi-shot EPI allows us mitigating  $T_2$  and  $T_2^*$  blurring and minimizing  $B_0$ -related distortion
- Combining the multiple shots can be challenging due to shot-to-shot phase variation.
- Advanced image reconstruction method is needed at high acceleration factor.


# Approach

- The recent zero-shot self-supervised learning (ZS-SSL) technique reconstructs images using scan-specific neural networks trained without additional training datasets.
- We employ ZS-SSL approach for joint reconstruction of accelerated multi-shot diffusion MRI.

# Approach



Unrolled network  
employed in our work



# Approach

- We train one network across all diffusion directions instead of training multiple networks individually.
- We use magnitude constraint (MC) of two shot assuming the two shots have similar signal intensity.

$$\underbrace{\left( \frac{\|u-v\|_1}{\|v\|_1} + \frac{\|u-v\|_2}{\|v\|_2} \right)}_{\text{k-space loss}} + \underbrace{\gamma \left| |s_1| - |s_2| \right|}_{\text{Magnitude constraint}}$$

k-space loss

Magnitude constraint

@  $R = 5$ ,  $b = 1k$ , 32 directions,  $1 \times 1 \times 4 \text{mm}^3$


SENSE

Single shot  
ZS-SSL

2 shot  
LORAKS

Proposed

1<sup>st</sup> shot




30.35 %

16.53 %

15.12 %

11.82 %

2<sup>nd</sup> shot



NRMSE:

24.60 %

15.53 %

14.41 %

13.42 %

# Summary

- We introduced a ZS-SSL method for multishot diffusion MRI reconstruction.
- Our method can yield better reconstruction than state-of-the-art LORAKS reconstruction.
- Proposed model can be trained for all directions at once, with an additional benefit of reduced training time. (2.5-fold)
- The magnitude constraint improves reconstruction performance.